



# DELIVERS

The One-Step Pregnancy Test  
Suitable for Urine & Serum

**ITS SENSITIVE ...**  
25 mIU/L

**ITS USER FRIENDLY ...**  
Just ADD 6 drops of either URINE or SERUM and wait for results to develop.

**ITS SPECIFIC ...**  
Other hormones and drugs tested do not present any interference on the test result.

**ITS STABLE ...**  
Over one year shelf life at room temperature storage.

**ITS ACCURATE ...**  
The internal control assures that a negative result represents a true negative test. The coloured lines in windows B & C allow easy reading on positive results.

# ENTREGA

El Examen de Embarazo de un Solo Paso  
Adecuado Para Orina Y Serum

**ES SENSIBLE ...**  
25 mIU/L

**ES FACIL DE USAR ...**  
Solamente AÑADA 6 gotas de ORINA o de SERU y espere por los resultados.

**ES ESPECIFICO ...**  
Otras hormonas y drogas examinadas no presentan interferencia alguna en el resultado del examen.

**ES ESTABLE ...**  
Más de un año de duración en el anaquel almacenado a temperatura ambiente.

**ES EXACTO ...**  
El control interno asegura que un resultado negativo representa un verdadero examen negativo. Las líneas de color en las ventanas B & C permiten una fácil lectura en los resultados positivos.

## PREG SLIDE (URINE)

DIRECT, LATEX M

CAT. NO.	PACKAGING	NO. TESTS
9000-1	LIQUID	100
9000-2	LIQUID	50

CONTROL INCLUDED

### TECHNICAL INFORMATION

STORAGE TEMP: 2 - 8° C  
SAMPLE: URINE  
SENSITIVITY: > 300 mIU/mL (0.3 IUHCG)

## PREG CUBE (SERUM / URINE)

DIRECT, CASSETTE M

CAT. NO.	PACKAGING	NO. TESTS
9100-1	CASSETTE	5
9100-2	CASSETTE	25

CONTROL INCLUDED

### TECHNICAL INFORMATION

STORAGE TEMP: 2 - 8° C  
SAMPLE: SERUM / URINE  
SENSITIVITY: > 300 mIU/mL (0.3 IUHCG)

## PREG STRIP (SERUM / URINE)

DIRECT, STRIP M

CAT. NO.	PACKAGING	NO. TESTS
9200-1	STRIPS	25

CONTROL INCLUDED

### TECHNICAL INFORMATION

STORAGE TEMP: 2 - 8° C  
SAMPLE: SERUM / URINE  
SENSITIVITY: > 40 mIU/mL