


Specialty Quick Tests

The specialty quick tests are offered in three formats: one format is a plastic cassette, with a membrane strip exposed; the second is a membrane strip and the third is a laminated, paper test card. In the first two formats, the sample migrates chromatographically on the membrane by capillary action. For a positive result, a color band with the specific conjugate complex will form on the membrane. Absence of this colored band in the test region suggests a negative result. To serve as a procedural control, a colored band at control region always appears in the test area to show the test works properly. All cassette cards and strips are individually packaged in foil. In the third format the tests are preformed using a laminated, paper test card and "latex" technology.



COCAINE

TECHNICAL INFORMATION

CAT. NO.	PACKAGING	NO. TESTS
S7765	STRIPS	25



STORAGE TEMP: 4 - 30° C
 SAMPLE: URINE
 SENSITIVITY: >300 ug/mL

THC

TECHNICAL INFORMATION

CAT. NO.	PACKAGING	NO. TESTS
S7770	STRIPS	25



STORAGE TEMP: 4 - 30° C
 SAMPLE: URINE

HBsAg

TECHNICAL INFORMATION

CAT. NO.	PACKAGING	NO. TESTS
S7706	CASSETTE	25



STORAGE TEMP: 2 - 8° C
 SAMPLE: SERUM

NITROPRUSSIDE SET

TECHNICAL INFORMATION

CAT. NO.	PACKAGING	NO. TESTS
N-110	LIQUID	500



STORAGE TEMP: 2 - 8° C
 SAMPLE: SERUM
 SENSITIVITY: >50 ng/dL

HIV

TECHNICAL INFORMATION

CAT. NO.	PACKAGING	NO. TESTS
S9000-1	CASSETTE	10



STORAGE TEMP: 2 - 8° C
 SAMPLE: SERUM, PLASMA
 OR WHOLE BLOOD

SYPHILIS (RPR)

TECHNICAL INFORMATION

CAT. NO.	PACKAGING	NO. TESTS
S8000-1	CARDS	200
S8000-2	CARDS	500

3 LEVEL
 CONTROL
 INCLUDED

STORAGE TEMP: 2 - 8° C
 SAMPLE: SERUM OR PLASMA